

Smart Oxygen Backpack

Smart Oxygen Backpack

The inside configuration features two inch nylon straps to secure a "D" or jumbo "D" size oxygen tank. A large utility pouch features a clear vinyl window to view contents. Includes quick-draw bagvalve mask mesh pouch. Both pouches are removable. The inside of the front panel features various pockets, elastic loops, and a protective zippered vinyl cover to keep supplies organized and clean. All packed 4/case. Size: 20.5" x 16" x 5"

ITEM #	DESCRIPTION	QTY	PRICE
12019A-G	Polyester, Green	Each	
12019A-N	Polyester, Navy	Each	
12019A-OR	Polyester, Orange	Each	
12019A-R	Polyester, Red	Each	

ITEM #	DESCRIPTION	QTY	PRICE
12019C-RB	Nylon, Royal	Each	
12019C-R	Nylon, Red	Each	

