

Heavy Duty Body Bag with Curved Zipper

www.basic-medicalsupply.com

Heavy Duty Body Bag with Curved Zipper

This body bag designed by Basic Medical Supply. Body bag has 6 strong handles, 3 each side. This is commonly referred to as a Human Remains Pouch, Human Remains Bag, Human Remains Pouch or a Body Bag. Material: 18 mil HD vinyl coated scrim. Static Lift Tested to 550 LB Capacity. Each Bag contains 3 Toe Tags. All Packed 10/case.

ITEM #	DESCRIPTION	QTY	PRICE
17011	Complete Impervious Fabric, Black, 36" x 96"	Each	
17012	Complete Impervious Fabric, Black, 36" x 96"	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com