

## I.V. Stand,w/5 wheels

www.basic-medicalsupply.com

### I.V. Stand,w/5 wheels

This economy IV Stand is constructed of a chrome-plated steel post and an epoxy coated steel stand. The knobless extension is secured by a durable friction wheel and extends from 46 to 84 inches. The twin-wheeled caster configuration ensures fluid, easy movement.

ITEM #	DESCRIPTION	QTY	PRICE
MS-IV200	I. V. Stand,w/5 wheels_ Stainless steel	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com