

I.V. Arm Board

I.V. Arm Board

Disposable I.V. Armboards are made with a strong, corrugated core and thick, comfortable foam padding. Covered in vinyl with a washable satin inish, these boards are dielectrically sealed to prevent bacterial luid iniltration.

ITEM #	DESCRIPTION	QTY	PRICE
MS-IVAB318D	I. V. Arm Board_Color:White_Size:3"x18"	Each	
MS-IVAB315D	I. V. Arm Board_Color:White_Size:3"x15"	Each	
MS-IVAB312D	I. V. Arm Board_Color:White_Size:3"x12"	Each	
MS-IVAB39D	I. V. Arm Board_Color:White_Size:3"x9"	Each	

