

Standard Stair Chair(301)

www.basic-medicalsupply.com

Standard Stair Chair(301)

The Basic Stair Chair with Vinyl Cover makes maneuvering in cramped quarters easier for EMS professionals. This chair features 2-5" rear wheels that make for easy transport on most surfaces, as well as front and back handles for easier handling. Patients are secured during transport with 3 quick-release restraints. Includes a durable, easy-to-clean plastic back and seat, and vinyl cover. Load Capacity: 350 lbs, Size: 36" x 21" x 28", Weight: 19 lbs, Folded: 8".

ITEM #	DESCRIPTION	QTY	PRICE
MS-STA301	Orange, Aluminum	Each	
MS-STA302	Blue, 650X520X900mm, Aluminum	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com