

## Standard Stair Chair(321)

www.basic-medicalsupply.com

### Standard Stair Chair(321)

The stretcher is mainly used for high building's lift carrying patients and wounded person. The structure can be folded and the stretcher structure could get smaller. It has four wheels. The back of the stretcher is designed with four folded handles. Under the front of the stretcher with two flexible uplift rods. The stretcher has pieces of belts to approve the patient's safety during transport process. It is made of high-strength aluminum alloy materials. It is characterized by its being light-weighted, small-sized, easy-carried, use-safely and easy for sterilization and clean.

Product Size(L\*W\*H):700x520x910mm

Carton Size:940\*580\*180mm

Load Limit:159kg

N.W.:13kg,G.W.:8kG

ITEM #	DESCRIPTION	QTY	PRICE
MS-STA321	Black, Aluminum	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com