

Standard Stair Chair(701)

Standard Stair Chair(701)

The Basic Stair Chair with Vinyl Cover makes maneuvering in cramped quarters easier for EMS professionals. his chair features 2-5" rear wheels that make for easy transport on most surfaces, as well as front and back handles for easier handling. Patients are secured during transport with 3 quick-release restraints. Includes a durable, easy-to-clean plastic back and seat, and vinyl cover. Load Capacity:350 lbs, Size: 36" x 21" x 28", Weight: 19 lbs, Folded: 8".

ITEM #	DESCRIPTION	QTY	PRICE
MS-STA701	Black, Size:900X460X480mm, Material: Aluminum, 1pcs/Case	Each	

