

Plastic Scoop Stretcher

www.basic-medicalsupply.com

Plastic Scoop Stretcher

Designed to gently maneuver stretcher under patient without rolling or lifting. Immobilization of the patient in the position found, minimizes the risk of complicating the existing injuries. The center of the stretcher can be opened to allow the patient to be X-rayed while secured on the stretcher. Folds for easy storing and separates in half during application and removal. Constructed of impact-resistant, light-weight and easy-to-clean polyethylene plastic. The head section is recessed to provide proper cervical alignment, allowing the patient's spine to remain in the position found. Features: X-Ray Translucent, Lightweight Plastic, Included: (3) Patient Restraint Straps, Adjustable Length, Dimensions: 80" L x 16.5" W x 2.75" D, Load Capacity 400 lbs. Weight: 17 lbs.

ITEM #	DESCRIPTION	QTY	PRICE
MS-SCP124	Yellow, 1pcs/Case	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com