

Folding Stretcher, Aluminum

www.basic-medicalsupply.com

Folding Stretcher, Aluminum

Folding stretcher strengthened by two steel bars that allow its use even for heavy transport. It is foldaway for an easy location. This stretcher is intensified stretcher, with flexible handle and outrigger, can be folded by lengthwise. It characterized by its being easy-carried, use-safely and easy for sterilization and clean.

Product Size (L*W*H):1850*480*30mm

Folded Size (L*W*H):460*480*130mm

Carton Size (2pcs): 960*530*160mm

Load Limit: 159kg

N.W.:4k gG.W.: 9kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL007	Size: 1850*480*30mm, Material: Aluminum	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com