

Folding Stretcher, Aluminum

www.basic-medicalsupply.com

Folding Stretcher, Aluminum

This separable first-aid ambulance stretcher is mainly used for transporting patients with fractures. Both ends of the stretcher have hinge-type clutches to separate the stretcher into a right and left part. This allows for the patient to be quickly and gently placed onto the stretcher, operation table or bed. The stretcher is then withdrawn from beneath the patient.

Product Size (L×W×H):185×48×20cm

Folding Size (L×W×H):46×48×13cm

N.W.:7 kg

Load Bearing: less than 159kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL008	Size: 185×48×20cm, Material: Aluminum	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com