

## Easy Fold Wheeled Stretcher,w/Adjustable Back Rest

### Easy Fold Wheeled Stretcher,w/Adjustable Back Rest

Designed to gently maneuver stretcher under patient without rolling or lifting. Immobilization of the patient in the position found, minimizes the risk of complicating the existing injuries. The center of the stretcher can be opened to allow the patient to be X-rayed while secured on the stretcher. Folds for easy storing and separates in half during application and removal. Features: Lightweight Aluminum, Included 3 Patient Restraint Straps, Adjustable Length, Dimensions: 66 1/4" L x 17 1/2" W x 2 5/8" D, Dimensions Folded: 49 1/2" L x 3 1/2" D x Adjustable to 80", Load Capacity 450 lbs.

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL010	Green, Size: 66-1/4"Lx17-1/2"Wx2-5/8", 1pcs/Case	Each	
MS-FL110	Orange, Size: 66-1/4"Lx17-1/2"Wx2-5/8", 1pcs/Case	Each	

