

Basket Stretcher(3000)

Basket Stretcher(3000)

The stretcher is made with an eye on a very wide range of emergency situations, such as in the mountains, air and at sea. It is sturdy and flexible to use, it is quick reliable fittings enable first-aid personnel to operate rapidly and safely. With its special sling equipment, the stretcher is ideal for lifting and transport by helicopter. It has an adjustable feet-securing mechanism, safety belts and a mattress. All the materials used are non-toxic, non-polluting, fire-proof, wear-resistant and anti-corrosion. The stretcher has 4 wheels which can move on the ground, reducing the effort of the operator.

Product Size (L*W*H): 2180*640*300mm

Carton Size: 2220*680*330mm

Load Limit: 272kg N.W.:11kg G.W.:20kg

Item # MS-BS3000

ITEM #	DESCRIPTION	QTY	PRICE
MS-BS3000	Orange, Size: 2180*640*300mm	Each	

