

Basket Stretcher(4000)

Basket Stretcher(4000)

The stretcher is made with an eye on a very wide range of special emergency situations ,such as in the mountains,air and at sea. It is sturdy and flexible to use, its quick reliable fittings enable first-aid person to operate it rapidly and safely .With its special sling equipment, the stretcher is ideal for lifted and transported by helicopter. It has an adjustable feet-securing mechanism, two safety belts and a mattress.All the materials used are mpm-toxic,non-pulluting,fire-proof.Weafresistant and anticorrosion The stretcher can be float on the sea ,buoyancy,120kg.

Product Size(L*W*H)

2180*1090*200mm

Carton Size:2200*680*230mm

Load Limit:272kg

N.W.:11kg G.W.:20kg

Item#MS-BS4000

ITEM #	DESCRIPTION	QTY	PRICE
MS-BS4000	Orange, Size: 2180*1090*200	Each	

