

## Basket Stretcher(5000)

### Basket Stretcher(5000)

The stretcher is made with an eye on a very wide range of special emergency situations ,such as in the mountains,air and at sea.It is sturdy and flexible to use, its quick reliable fittings enable first-aid person to operate it rapidly and safely. With its special sling equipment, the stretcher is ideal for lifted and transported by helicopter. It has an adjustable feet-securing mechanism,two safety belts and a mattress. The basket stretcher is made of stainless steel tube and net, it will not rust and potentially cause additional harm to the patient.It is stronger than common basket stretcher.It can be used together with the Spine Boards.

Product Size (L\*W\*H): 2200\*620\*190mm

Carton Size: 2220\*650\*210mm

Load Limit: 408kg N.W.:11kg G.W.:17kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-BS5000	Size: 2200*620*190mm	Each	

