

## Automatic Loading Stretcher(003)

### Automatic Loading Stretcher(003)

It is made of high-strength aluminum alloy; it is characterized by its being light-weighted, use-safely and easy for sterilization and cleaning. The stretcher has four swivel wheels. It's mainly used in hospitals, ambulance vehicles for carrying patients and wounded person Product Size (L×W×H)

1950x550x108mm □ 340mm □

Maximum Angle of the Back 60°

Carton Size 200×620×380mm

Load Limit 250kg

N.W 50kg G.W 55kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-ASL003	Orange, Size:1950x550x108mm (340mm) , High-strength Aluminum, 1pcs/Case	Each	

