

Automatic Loading Stretcher(005)

Automatic Loading Stretcher(005)

It is made of high-strength aluminum alloy; it is characterized by its being light-weighted, use-safely and easy for sterilization and cleaning. The stretcher has four swivel wheels. It's mainly used in hospitals, ambulance vehicles for carrying patients and wounded person Product Size (L×W×H) 1930x550x830mm(260mm)

Maximum Angle of the Back 60°

Carton Size 198×620×290mm

Load Limit 180kg

N.W 38kg G.W 43kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-ASL005	Orange, 1930x550x830mm(260mm), High-strength Aluminum, 1pcs/Case	Each	

