

Automatic Loading Stretcher(006)

Automatic Loading Stretcher(006)

It is made of high-strength aluminum alloy; it is characterized by its being light-weighted, use-safely and easy for sterilization and cleaning. The stretcher has four swivel wheels. It's mainly used in hospitals, ambulance vehicles for carrying patients and wounded person Product Size (L×W×H) 1980x550x112mm(280mm)

Maximum Angle of the Back 60°

Carton Size 205×620×290mm

Load Limit 250kg

N.W 50kg G.W 55kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-ASL006	Orange, 1980x550x112mm(280mm), High-strength Aluminum, 1pcs/Case	Each	

