

Half-body CPR Manikin

Half-body CPR Manikin

Features :

1. Simulating standard open airway
2. External breast compression: display device and alarm device
 - a. indicator light display of correct and wrong compression ; alarm of wrong compression
 - b. intensity display of correct (4-5 cm) and wrong (<4cm or >5cm) compression; alarm of wrong compression .
3. Artificial respiration (inhalation) : display device and alarm device
 - a. Inhalation <500-600ml or > 600ml ,wrong indicator light display and alarm prompting ;Inhalation between 500-600ml right indicator light display
 - b. Indicator light display open airway
 - c. Inhalation too quickly or too much result in air entering into stomach ; indicator light display and alarm prompting .
- 4 Ratio of compression and artificial respiration: 30:2 (one or two person).
5. Operating cycle: one cycle includes five times of 30:2 ratio of compression and artificial respiration.
6. Operation frequency: 100 times per minute (the latest international standard)
7. Operation methods: exercise operation
8. Examination of pupil response: mydriasis and myosis
9. Examination of carotid response : pinch the pressure ball by hand and simulate carotid pulse
10. Working Conditions: output power is 6V transformed from 220V external power through manostat.

ITEM #	DESCRIPTION	QTY	PRICE
MS-TM200	Half-body CPR Manikin_Material:PVC	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com